

Acuven™

Rapid COVID-19 Testing with confidence

Acuven™ Rapid COVID-19
Testing with confidence



Acuven™

COVID-19 ready

Rapid test solution

First line COVID19- testing

Reliable and precise

Back boning the health systems

Testing is critical to controlling the coronavirus and eventually easing restrictions that have halted daily life. Currently, almost all testing in hospitals, clinics and drive-thru sites uses the first line mass screening method using rapid test cassettes, to help doctors detect and treat people with active COVID-19.

- Rapid-screening within 10-15 minutes.
- High detection efficiency: simultaneous monitoring of IgM and IgG.
- Detection without any testing equipment.
- Easy to operate, and is compatible with serum/ whole blood/ plasma.
- Room-temperature storage.

STERILE A

Sterilized using aseptic processing technique



Caution, consult accompanying documents



Do not reuse

IVD



Meroven LLC USA - 7901 4th St. N STE 300, St. Petersburg, FL., USA 33702

www.meroven.us

E-mail: info@meroven.com